

MICHAEL A. TRUJILLO

michael.anthony.trujillo.jr@gmail.com
720.620.7438

REFERENCES

Marty Chenery 303.987.9797

Joel Garret 303.746.5337

Kai Metcalfe 720.317.3637

Tom Kidwell 317.459.5816

Ryan Rizzio 847.951.2853

EXPERIENCE

Manufacturing Associate II

Quicksilver Scientific / Nov 2020 – Feb 2023

- Responsible for handling, weighing, and processing wet and dry ingredients
- Maintained the cleanliness and general upkeep of the controlled work area
- Reviewed and responded to batch records
- Utilized fast problem solving to maintain manufacturing consistency, efficiency, and operational standards
- Responsible for breaking down line and equipment and cleaning and sanitizing
- Operated homogenizers, processed vitamin supplements and CBD products
- Experienced with particle sizing software
- Capable of operation of large pressurized vessels

Bottling Technician (second shift)

Quicksilver Scientific / Mar 2020 – Nov 2020

- Set up fill line according to product in use, bottle, cap size
- Broke down line and equipment and conducted cleaning and sanitization
- Maintained CGMP documentation in logbooks
- Performed preventative maintenance and troubleshooting for equipment per SOP
- Responsible for prepping bottling and packaging equipment for next day production goals
- Followed pass down for setting up bottling equipment for 1st shift production

CBD Encapsulation Lead

Functional Remedies / Dec 2017 – Oct 2019

- Trained over ten technicians in two years and led the encapsulation department team which consisted of five technicians
- Conducted CBD oil extractions, including homogenization and dilutions, to obtain varying product potencies as well as operated various machinery for manufacturing CBD products such as capsules, tinctures, and salves
- Conducted calibration and maintenance of peristaltic pumps and pneumatic pistons for tinctures and capsules, respectively.
- CGMP experience with four years of aseptic technique practice.

Quality Inspector

Thomas & Betts / Jan 2017 – Apr 2017

- Conducted quality inspections of plant manufactured electrical components, including transformer and junction box elbows bushings and coupling devices.
- Used visual inspection techniques, gauges, and x-ray inspections to pass or fail parts based on measurement, compound, and trimming tolerances according to SOP.

Wash and Prep Technician

Oso Bio Pharmaceuticals / Aug 2014 – May 2016

- Responsible for sterilizing and prepping equipment for clean room use, loading and unloading autoclave ovens while ensuring proper operations in an A/B space clean room environment.